

BOGHI Engineering & Software Technologies Private Limited

Branch Off: 4th Floor; CoKarma, Veer Chambers, Begumpet, Secunderabad, Telangana State-500016. Cell:9121003011

Regd. Off: Plot No.11, Anand Nagar-1, Bandlaguda, Nagole, Hyderabad, Telangana State-500068. Cell: 9347006933/9704010147

www.boghiengineering.com



BLACK & WHITE VISUAL INSPECTION BACKGROUND WITH LIGHT & MAGNIFYING GLASS (AS Make-USA) **(WALL HANGING TYPE)**



A unique product within the cleanroom industry - the Black & White Visual Inspection Background to provide a solution to all your inspection needs - one-of-a-kind inspection backgrounds are specifically designed and manufactured with the USP <790> and USP <797> visual inspection requirements in mind. This background is made to be attached to a wall, allowing for a consistent inspection setup & designed to make foreign particles stand out clearly in sterile containers and solutions during a manual visual inspection. With an attached LED lightbulb & adjustable magnifying glass, the background can be used in nearly any type of lighting and still be effective. Packaging for Visual Inspection Background is sold individually (includes wall-mount & screws for easy installation).

Specifications of Visual Inspection Background:

- **Material:** Made of a chemically resistant material for easy cleaning/maintenance
- **Chemical Resistance:** Yes
- **Style:** Wall-hanging
- **Size:** Standard
- **Dimensions:** 15.5" H x 16" W

Lighting Specifications:

- **Light:** Standard (LED)
- **Voltage:** 120 V (US Standard)
- **Lumen:** 360 Lumens
- **Brightness (Lux):** 2000 lux
- **Colour Rating (Kelvin):** 4,200 K
- **Rated Life:** 20,000 Hours
- **Dimensions of light fixture:** 13.5" L x 1.7" W

Compliances:

- Each inspection background is specifically designed and manufactured with the USP <790> and USP <797> visual inspection requirements in mind.

Applications:

- Designed to make foreign particles stand out clearly in sterile containers and solutions.

BOGHI Engineering & Software Technologies Private Limited

Branch Off: 4th Floor; CoKarma, Veer Chambers, Begumpet, Secunderabad, Telangana State-500016. Cell:9121003011

Regd. Off: Plot No.11, Anand Nagar-1, Bandlaguda, Nagole, Hyderabad, Telangana State-500068. Cell: 9347006933/9704010147

www.boghiengineering.com



BLACK & WHITE VISUAL INSPECTION BACKGROUND WITH LIGHT (AS Make-USA) (TABLE TOP TYPE)



Specifications of Visual Inspection Background:

- **Material:** Made of a chemically resistant material for easy cleaning/maintenance
- **Chemical Resistance:** Yes
- **Style:** Tabletop
- **Size:** Standard
- **Dimensions:** 15.5" H x 16" W

Lighting Specifications:

- **Light:** Standard (LED)
- **Voltage:** 120 V (US Standard)
- **Lumen:** 360 lumens
- **Brightness (Lux):** 2000 lux
- **Colour Rating (Kelvin):** 4,200 K
- **Rated Life:** 20,000 Hours
- **Dimensions of light fixture:** 13.5" L x 1.7" W

Compliances:

- Each inspection background is specifically designed and manufactured with the USP <790> and USP <797> visual inspection requirements in mind.

Applications:

- Designed to make foreign particles stand out clearly in sterile containers and solutions.

BOGHI Engineering & Software Technologies Private Limited

Branch Off: 4th Floor; CoKarma, Veer Chambers, Begumpet, Secunderabad, Telangana State-500016. Cell:9121003011

Regd. Off: Plot No.11, Anand Nagar-1, Bandlaguda, Nagole, Hyderabad, Telangana State-500068. Cell: 9347006933/9704010147

www.boghiengineering.com



BLACK & WHITE VISUAL INSPECTION BACKGROUND (AS Make-USA) (TABLETOP TYPE)



Specifications of Visual Inspection Background:

- **Material:** Made of a chemically resistant material for easy cleaning/maintenance
- **Chemical Resistance:** Yes
- **Style:** Tabletop
- **Size:** Standard
- **Dimensions:** 15.5" H x 16" W

Lighting Specifications:

- No Lighting included

Compliances:

- Each inspection background is specifically designed and manufactured with the USP <790> and USP <797> visual inspection requirements in mind.

Applications:

- Designed to make foreign particles stand out clearly in sterile containers and solutions.